



Practicon® Floss-Hands™

Floss Aids

Instructions

Brushing alone cleans only 60% of tooth surfaces. Flossing removes food particles and plaque from between teeth and under the gumline where a toothbrush can't reach.


Plaque buildup leads to bad breath, bleeding gums, periodontal disease, and tooth decay.

Securing floss:

1. Place one end of a 12" length of dental floss across button slot and wind floss around the button.
2. Pull floss tight across the floss span opening, following the grooves.
3. Pull floss tight across the button slot and wind around the button again several times.

Using Floss-Hands:

1. Gently glide (never snap) the floss between adjacent teeth and under the gumline as far as possible.
2. Holding floss against the tooth surface, sweep it away from the gums. Floss both sides of each tooth and behind back teeth.



Manufactured for

PRACTICON

PRACTICAL INNOVATIONS FOR DENTISTRY

1-800-959-9505 • Fax 1-800-278-0885

practicon.com • Greenville, NC 27834

Made in China • IFU 7039740_2