

Q77272  
Q28572  
U87295

ID# 4386

DATE: 1/97

MATERIAL SAFETY DATA SHEET

**PRODUCT IDENTIFICATION**

PRODUCT/TRADE NAME: SANI CLOTH GERMICIDAL TOWELETTES  
MANUFACTURER'S NAME: PDI - PROFESSIONAL DISPOSABLES, INC.  
THE HEALTHCARE DIVISION OF  
NICE-PAK PRODUCTS, INC.  
ADDRESS: 2 NICE-PAK PARK  
ORANGEBURG NY 10962-1376

**HAZARD RATING (NFPA)**

HEALTH: 0  
FLAMMABILITY: 2  
REACTIVITY: 0  
SPECIFIC: NONE

EMERGENCY OR INFORMATION TELEPHONE NO: 914-365-1700 (M-F DAYTIME) AT OTHER TIMES, CONTACT THE LOCAL POISON CONTROL CENTER

CHEMICAL NAME: PROPRIETARY MIXTURE

**I HAZARDOUS INGREDIENTS** PER 29 CFR 1910.1200

Hazardous Ingredients	%	ACGIH TLV	CAS NUMBER
ISOPROPYL ALCOHOL	10	400 ppm	67-63-0
QUATERNARY AMMONIUM CMPD	<1	N/A	27479-28-3

**III-PHYSICAL/CHEMICAL CHARACTERISTICS**

COLOR/ODOR/APPEARANCE: APPLICATOR SATURATED WITH CLEAR LIQUID WITH ALCOHOL ODOR  
BOILING POINT: 181° F  
FLASH POINT: 103° F  
VAPOR DENSITY N/A  
EVAPORATION RATE: N/A  
SOLUBILITY IN WATER: COMPLETE  
SPECIFIC GRAVITY (H<sub>2</sub>O = 1): .983

**IV FIRE & EXPLOSION HAZARD DATA**

FLASH POINT(Method Used): 103° F TAG CLOSED CUP : LEL: N/A UEL: N/A  
EXTINGUISHING MEDIA: WATER, FOAM

SPECIAL FIRE FIGHTING PROCEDURES: TREAT PRODUCT AS COMBUSTIBLE

UNUSUAL FIRE AND EXPLOSION HAZARDS: NONE

**V-REACTIVITY DATA**

STABILITY: STABLE  
CONDITIONS TO AVOID : NONE  
INCOMPATIBILITY: ANIONIC SURFACTANT, SOAPS  
HAZARDOUS DECOMPOSITION OR BYPRODUCTS: NONE KNOWN  
POLYMERIZATION: WILL NOT OCCUR  
CONDITIONS TO AVOID: NONE KNOWN  
PRODUCT DOES NOT REACT WITH FORMALDEHYDE

SANI CLOTH

**VI-HEALTH HAZARD DATA**

EFFECTS OF OVER EXPOSURE:

SKIN : MAY CAUSE SKIN IRRITATION  
EYES : MAY CAUSE IRRITATION  
INHALATION: NOT NORMAL ROUTE OF ENTRY  
INGESTION: NOT NORMAL ROUTE OF ENTRY

EMERGENCY AND FIRST AID PROCEDURES:

SKIN CONTACT: IF RASH OR IRRITATION DEVELOPS, DISCONTINUE USE  
EYE CONTACT: FLUSH WITH COOL WATER FOR 15 MINUTES  
INHALATION: REMOVE TO FRESH AIR  
INGESTION: CONSULT PHYSICIAN  
TARGET ORGANS; NONE KNOWN  
KNOWN CARCINOGEN: NONE KNOWN

**VII-SPILL AND DISPOSAL PROCEDURE**

SPILL CONTROL: PICK UP WITH ABSORBENT

WASTE DISPOSAL METHOD: TREAT AS FLAMMABLE LIQUID. OBEY LOCAL, STATE AND FEDERAL REGULATIONS

HANDLING AND STORAGE: AVOID HIGH HEAT AND SOURCES OF IGNITION

**VIII-CONTROL MEASURES/PROTECTION**

RESPIRATION: USE NIOSH APPROVED RESPIRATOR IF EXPOSURE LEVEL IS EXCEEDED  
VENTILATION: RECOMMENDED  
PROTECTIVE GLOVES: RECOMMENDED  
EYE PROTECTION: IF SPLASHING MAY OCCUR  
HYGIENIC PRACTICES: FOLLOW GOOD HOUSEKEEPING PRACTICES  
OTHER: NONE

**IX TRANSPORT/SHIPPING**

DOT SHIPPING NAME:  
TECHNICAL SHIPPING NAME:  
DOT SHIPPING CLASSIFICATION:  
DOT ID NO.:  
DOT LABEL; REQUIREMENTS:  
UN/NA NUMBER REGULATIONS:  
REPORTABLE QUANTITY:

**X DISCLAIMER**

THE INFORMATION FURNISHED HEREIN IS BELIEVED TO BE ACCURATE AND REPRESENTS THE BEST DATA CURRENTLY AVAILABLE TO US. NO WARRANTY, EXPRESSED OR IMPLIED IS MADE AND NICE-PAK PRODUCTS, INC. ASSUMES NO LEGAL RESPONSIBILITY OR LIABILITY RESULTING FROM ITS USE.